

INSTRUCTION SHEET SUBMISSION OF PATIENT BLOOD SAMPLE

- For red blood cell request: draw a minimum of TWO: Full 6 mL EDTA (purple top or pink top) tubes within 3 days of intended transfusion date.
- See table below for the sample requirements.

Transfusion History	# of Tubes to be drawn	Affix BBID to
New patient with no history of patient's blood type	 2 tubes drawn at different times by two separate individuals if possible Verification of patient identification must be made prior to collection of each sample 	Both tubes
Patient's historical blood type is available per LifeStream record	2 tubes can be drawn at the same time	Both tubes

- Ideally, the samples will be drawn by 2 people drawing 2 separate samples. If staffing is prohibitive, one staff member may draw two separate tubes.
- Affix a BBID number matching the one attached to the patient's wristband to the Request for Transfusion of Blood Components form.
- Apply printed label if available to each tube of blood with the following information:
 - Patient's full name, last name first, including middle initial, if any.
 - Patient Medical Record Number.
 - Date of birth.
 - Date and time of collection.
 - Initials of person drawing the sample.
- Verify patient's name, patient Medical Record Number, date of birth, and BBID prior to collection of each sample.
- Verify information on the Request for Transfusion of Blood Components Form with information on the tube label prior to submitting to LifeStream.
- Place all tubes in a plastic biohazard bag after collection. Attach the Request for Transfusion of Blood Components to the outside of the bag.
- Keep samples at room temperature or refrigerate samples between 1-10°C until ready to transport to LifeStream.
- Call LifeStream Reference Laboratory at (909) 386-6858 to request the pickup.

NOTE: If additional units of blood are needed for a patient and the previous specimen was drawn within 3 days of the next transfusion date, new samples are not required. Submit only a new request by fax.